

# VITAMIN B6

## Orderable – VITB6

Turnaround Time: 10 days

### Alternate Name(s):

Pyridoxal Phosphate  
Pyridoxal-5-Phosphate

### Specimen:

4.5 mL Green Lithium Heparin top or 4 mL K<sub>2</sub> or K<sub>3</sub> EDTA Lavender top Vacutainer tube

### Collection Information:

Fasting sample is recommended.

Patient must avoid dietary Vitamin B6 on day of collection.

Minimum volume required is 1 mL.

Protect sample from light by wrapping in foil.

### Reference Ranges:

20-96 nmol/L

### Special Processing:

Protect sample from light.

### Storage and Shipment:

Separate and freeze.



**Laboratory:**  
Toxicology/Special  
Chemistry Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
LC-MS/MS



**Test Schedule:**  
Monday--Friday  
0800-1600