

SCHISTOSOMA SEROLOGY

Orderable – PARASRO

Turnaround Time: Up to 10 days

Specimen Type and Collection Information:

Please review the Public Health Ontario Laboratory Test Information Index:

Schistosoma Serology

Specimen Type	Minimum Volume	Collection Kit
Whole blood or serum	5.0 mL whole blood or 1.0 mL serum	Gold 5mL Vacutainer (SST)



Laboratory:
Microbiology



Requisition:
[PUBLIC HEALTH
LABORATORY TEST
REQUISITION](#)

This requisition **MUST**
accompany the specimen



Method of Analysis:
Enzyme Linked
Immunosorbant Assay
(ELISA)



Test Schedule:
Referred out weekdays



Referred to:
Public Health Laboratory
1200 Commissioners Rd E
London, Ontario
519-455-9310

Critical Information Required:

Ensure the Public Health Laboratory Requisition includes the following:

- Onset date
- Symptoms

Specify **Schistosoma** in the **Specimen Description** field when placing order in PowerChart and on Public Health Laboratory Requisition.

Storage and Shipment:

The specimen must be refrigerated.
The specimen must be received within 24 hours of collection.