

APO
LIPOPROTEIN
B

Orderable - APOB

Turnaround Time: 14 days

Alternate Name(s):

Apo B
ApoB

Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer tube	2-10 years: 3 mL Red top Vacutainer tube
Lavender top (EDTA) tubes are also acceptable	
Light Green (Li-Heparin) top tubes are NOT acceptable	

Collection Information:

Minimum volume of serum required is 1 mL for adult or pediatric samples.

Reference Ranges:

0.59-1.46 g/L

Interpretive Comments:

Please note change in method as of Aug 5, 2021.

When triglycerides are >10 mmol/L, analysis is not accurate and will not be performed.

Storage and Shipment:

Store and ship serum frozen.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Nephelometry



Test Schedule:
Monday-Thursday



Referred Out Location:
[In-Common Laboratories](#)