

ALBUMIN, PLASMA

Orderable - ALB

Turn Around Time: 4 hours

STAT: 1 hour

Specimen:

Adult	Pediatric
4.5 mL Green (Lithium Heparin) top Vacutainer	0-2 years: 0.5 mL Green top Microtainer 2-10 years: 2 mL Green top Vacutainer tube

Collection Information:

Collect blood aseptically in a Vacutainer tube.

Minimum volume required for a pediatric sample is 0.2 mL

Reference Ranges:

Age	Range
0-4 days:	28-44 g/L
5 days-14 years:	38-54 g/L
15 years-18 years:	32-45 g/L
>18 years:	35-52 g/L



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Colorimetric endpoint
(BCG)



Test Schedule:
As required