

HEPATITIS B CORE ANTIBODY

Orderable - HpCoreAb

Turn Around Time: 48 hours

Alternate Name(s):

CORE
Anti HBc (IgG + IgM)

Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer tube	0-2 years: 0.5 mL Gold top Microtainer 2-10 years: 2 mL Gold top tube

Collection Information:

Collect blood aseptically in a Vacutainer tube.

Reference Ranges:

See report

Interpretive Comments:

Testing detects both IgG and IgM antibodies (total). This test can be a life-long Hepatitis marker that represents past infection with Hepatitis B.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Electrochemiluminescence
Immunoassay



Test Schedule:
24/7