



HEPATITIS B DIAGNOSTIC PANEL

Orderable - HPBPANEL

Turn Around Time: 48 hours

Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer tube	0-2 years: 0.5 mL Gold top Microtainer 2-10 years: 2 mL Gold top tube



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Electrochemiluminescence
Immunoassay



Test Schedule:
24/7

Collection Information:

Collect blood aseptically in a Vacutainer tube.

Interpretive Comments:

Panel includes Hepatitis B surface antigen, Hepatitis B surface antibody, Hepatitis B core antibody (IgG and IgM)