



HEPATITIS D ANTIBODY

Orderable - HpDAb

Turn Around Time: 30 days

Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer tube	0-2 years: 0.5 mL Gold top Microtainer 2-10 years: 2 mL Gold top tube



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION

Collection Information:

Collect blood aseptically in a Vacutainer tube.

Reference Ranges:

See report

Special Processing:

Specimen will only be processed if Hepatitis B surface antigen is positive.



Test Schedule:
Referred out Monday-
Friday between 0800-
1600 to Public Health
Laboratory



Referred Out Location:
[Public Health Laboratory](#)