

IODINE,  
URINE:  
Random or  
24-hour

## Orderable - IU/IU24

Turn Around Time: 10 days

### Specimen:

24-hour urine collected in an **unused** 24-hour urine container or random urine

### Collection Information:

A minimum volume of 15 mL of urine is required.

#### Collection of 24-Hour Urine Sample

Collect and transfer in metal-free container. Avoid mineral supplements for 5 days.

For random collections, indicate 'random'.

Please complete the laboratory label attached to the collection bottle to include patient's name, collection start time and date and collection end time and date.

### Reference Ranges:

<b>24-Hour urine:</b>
0.79-3.62 umol/d
Conversion factor ug/d x 0.0079

<b>Random urine:</b>
0.33-2.76 umol/L
conversion factor ug/dL x 0.0788



**Laboratory:**  
Core Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
ICP-TripleQ MS



**Test Schedule:**  
Referred out Monday-  
Thursday



**Referred Out Location:**  
[In-Common Laboratories](#)